
PCNP antibody

Product Information

Catalog No.:	FNab06218
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May be involved in cell cycle regulation.

Immunogen information


Immunogen:	PEST proteolytic signal containing nuclear protein
Synonyms:	PEST proteolytic signal-containing nuclear protein (PCNP, PEST-containing nuclear protein) PCNP
Observed MW:	24-28 kDa
Uniprot ID :	Q8WW12

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human lymphoma using FNab06218(PCNP antibody) at dilution of 1:50



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab06218(PCNP antibody) at dilution of 1:600