

PCNX2 antibody

Product Information

Catalog No.:	FNab06219
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	pecanex-like 2(Drosophila)
Synonyms:	Pecanex-like protein 2 Pecanex homolog protein 2 PCNX2 KIAA0435 PCNXL2
Observed MW:	236 kDa
Uniprot ID :	A6NKB5

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:1000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab06219(PCNXL2 antibody) at dilution of 1:500