

PCP4L1 antibody

Product Information

Catalog No.:	FNab06223
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	Purkinje cell protein 4 like 1
Synonyms:	Purkinje cell protein 4-like protein 1 (PCP4-like protein 1) PCP4L1 IQM1
Observed MW:	10 kDa
Uniprot ID :	A6NKN8

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded mouse brain tissue slide using FNab06223(PCP4L1 Antibody) at dilution of 1:400



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab06223(PCP4L1 Antibody) at dilution of 1:600