

## PCSK7 antibody

### Product Information

Catalog No.:	FNab06226
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Likely to represent a ubiquitous endoprotease activity within constitutive secretory pathways and capable of cleavage at the RXXX[KR]R consensus motif.

### Immunogen information


Immunogen:	proprotein convertase subtilisin/kexin type 7
Synonyms:	Proprotein convertase subtilisin/kexin type 7 Lymphoma proprotein convertase Prohormone convertase 7 Proprotein convertase 7 (PC7) Proprotein convertase 8 (PC8, hPC8) Subtilisin/kexin-like protease PC7 PCSK7 LPC PC7 PC8 SPC7
Observed MW:	68 kDa
Uniprot ID :	Q16549

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human colon cancer using FNab06226(PCSK7 antibody) at dilution of 1:50



mouse spleen tissue were subjected to SDS PAGE followed by western blot with FNab06226(PCSK7 Antibody) at dilution of 1:300