

## CDK16 antibody

### Product Information

Catalog No.:	FNab06228
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	PCTAIRE protein kinase 1
Synonyms:	Cyclin-dependent kinase 16 Cell division protein kinase 16 PCTAIRE-motif protein kinase 1 Serine/threonine-protein kinase PCTAIRE-1 CDK16 PCTAIRE1 PCTK1
Observed MW:	56-63 kDa
Uniprot ID :	Q00536

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:200-1:2000
Image:	



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab06228(PCTK1 antibody) at dilution of 1:600