

CDK18 antibody

Product Information

Catalog No.:	FNab06229
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	PCTAIRE protein kinase 3
Synonyms:	Cyclin-dependent kinase 18 Cell division protein kinase 18 PCTAIRE-motif protein kinase 3 Serine/threonine-protein kinase PCTAIRE-3 CDK18 PCTAIRE3 PCTK3
Observed MW:	54 kDa
Uniprot ID :	Q07002

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab06229(PCTAIRE3 antibody at dilution of 1:1000