

## PCTP antibody

### Product Information

Catalog No.:	FNab06230
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Catalyzes the transfer of phosphatidylcholine between membranes. Binds a single lipid molecule.

### Immunogen information


Immunogen:	phosphatidylcholine transfer protein
Synonyms:	Phosphatidylcholine transfer protein (PC-TP) START domain-containing protein 2 (StARD2) StAR-related lipid transfer protein 2 PCTP STARD2
Observed MW:	25-27 kDa
Uniprot ID :	Q9UKL6

### Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human ovary tumor using FNab06230(PCTP antibody) at dilution of 1:50



human liver tissue were subjected to SDS PAGE followed by western blot with FNab06230(PCTP antibody) at dilution of 1:300