

PCYOX1 antibody

Product Information

Catalog No.:	FNab06231
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Involved in the degradation of prenylated proteins. Cleaves the thioether bond of prenyl-L-cysteines, such as farnesylcysteine and geranylgeranylcysteine.

Immunogen information

Immunogen:	prenylcysteine oxidase 1
Synonyms:	Prenylcysteine oxidase 1 Prenylcysteine lyase PCYOX1 KIAA0908 PCL1
Observed MW:	None
Uniprot ID :	Q9UHG3

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human colon using FNab06231(PCYOX1 antibody) at dilution of 1:50