

PCYT1B antibody

Product Information

Catalog No.:	FNab06232
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Controls phosphatidylcholine synthesis.

Immunogen information

Immunogen:	phosphate cytidylyltransferase 1, choline, beta
Synonyms:	Choline-phosphate cytidylyltransferase B CCT-beta CTP:phosphocholine cytidylyltransferase B (CCT B, CT B) Phosphorylcholine transferase B PCYT1B CCTB
Observed MW:	42 kDa
Uniprot ID :	Q9Y5K3

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:1000
Image:	



human placenta tissue were subjected to SDS PAGE followed by western blot with FNab06232(PCYT1B antibody) at dilution of 1:500