

PD-1 antibody

Product Information

Catalog No.:	FNab06235
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	2E8
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	programmed cell death 1
Synonyms:	Programmed cell death protein 1 (Protein PD-1, hPD-1) PDCD1 PD1
Observed MW:	32 kDa, 50 kDa
Uniprot ID :	Q15116

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, FC
Recommended dilution:	WB: 1:2000-1:10000; IHC: 1:1000-1:5000
Image:	



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using FNab06235(PD-1/CD279 Antibody) at dilution of 1:1000



Raji cells were subjected to SDS PAGE followed by western blot with FNab06235(PD-1/CD279 Antibody) at dilution of 1:2000