

PDCD7 antibody

Product Information

Catalog No.:	FNab06244
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Promotes apoptosis when overexpressed.

Immunogen information

Immunogen:	programmed cell death 7
Synonyms:	Programmed cell death protein 7 ES18 (hES18) PDCD7
Observed MW:	55 kDa
Uniprot ID :	Q8N8D1

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab06244(PDCD7 antibody) at dilution of 1:3000