

PDCL3 antibody

Product Information

Catalog No.:	FNab06247
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Acts as a chaperone for the angiogenic VEGF receptor kDaR/VEGFR2, controlling its abundance and inhibiting its ubiquitination and degradation. Modulates the activation of caspases during apoptosis. Is a substrate for Orgyia pseudotsugata multicapsid polyhedrosis virus(OpMNPV) IAP-mediated ubiquitination.

Immunogen information

Immunogen:	phosducin-like 3
Synonyms:	Phosducin-like protein 3 HTPHLP PhPL3 Viral IAP-associated factor 1 (VIAF-1) PDCL3 PhLP2A VIAF1
Observed MW:	25 kDa
Uniprot ID :	Q9H2J4

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



A2780 cells were subjected to SDS PAGE followed by western blot with FNab06247(PDCL3 antibody) at dilution of 1:3000