

PDE4C antibody

Product Information

Catalog No.:	FNab06254
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Hydrolyzes the second messenger cAMP, which is a key regulator of many important physiological processes.

Immunogen information

Immunogen:	phosphodiesterase 4C, cAMP-specific(phosphodiesterase E1 dunce homolog, Drosophila)
Synonyms:	3',5'-cyclic-AMP phosphodiesterase 4C DPDE1 PDE21 cAMP-specific phosphodiesterase 4C PDE4C DPDE1
Observed MW:	80 kDa
Uniprot ID :	Q08493

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500 - 1:2000
Image:	



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab06254(PDE4C antibody) at dilution of 1:1000