

PDE6B antibody

Product Information

Catalog No.:	FNab06258
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background

This protein participates in processes of transmission and amplification of the visual signal. Necessary for the formation of a functional phosphodiesterase holoenzyme.

Immunogen information

Immunogen:	phosphodiesterase 6B, cGMP-specific, rod, beta
Synonyms:	Rod cGMP-specific 3',5'-cyclic phosphodiesterase subunit beta (GMP-PDE beta) PDE6B PDEB
Observed MW:	None
Uniprot ID :	P35913

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA