

PDE6B antibody

Product Information

Catalog No.:	FNab06259
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


This protein participates in processes of transmission and amplification of the visual signal. Necessary for the formation of a functional phosphodiesterase holoenzyme.

Immunogen information

Immunogen:	phosphodiesterase 6B, cGMP-specific, rod, beta
Synonyms:	Rod cGMP-specific 3',5'-cyclic phosphodiesterase subunit beta (GMP-PDE beta) PDE6B PDEB
Observed MW:	98 kDa
Uniprot ID :	P35913

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



A549 cells were subjected to SDS PAGE followed by western blot with FNab06259(PDE6B antibody) at dilution of 1:3000