

PDE6C antibody

Product Information

Catalog No.:	FNab06260
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	phosphodiesterase 6C, cGMP-specific, cone, alpha prime
Synonyms:	Cone cGMP-specific 3',5'-cyclic phosphodiesterase subunit alpha' cGMP phosphodiesterase 6C PDE6C PDEA2
Observed MW:	55 kDa, 99 kDa
Uniprot ID :	P51160

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



Y79 cells were subjected to SDS PAGE followed by western blot with FNab06260(PDE6C antibody) at dilution of 1:3000