

## TYMP antibody

### Product Information

Catalog No.:	FNab06264
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	thymidine phosphorylase
Synonyms:	Thymidine phosphorylase (TP) Gliostatin Platelet-derived endothelial cell growth factor (PD-ECGF) TdRPase TYMP ECGF1
Observed MW:	50 kDa
Uniprot ID :	P19971

### Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human ovary tumor using FNab06264(PD-ECGF antibody) at dilution of 1:50



human placenta tissue were subjected to SDS PAGE followed by western blot with FNab06264(PD-ECGF antibody) at dilution of 1:300