

PDLIM1 antibody

Product Information


Catalog No.:	FNab06282
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	PDZ and LIM domain 1
Synonyms:	PDZ and LIM domain protein 1 C-terminal LIM domain protein 1 Elfin LIM domain protein CLP-36 PDLIM1 CLIM1 CLP36
Observed MW:	50 kDa
Uniprot ID :	O00151

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, IF, WB, FC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human pancreas cancer using FNab06282(PDLIM1 antibody) at dilution of 1:50



mouse lung tissue were subjected to SDS PAGE followed by western blot with FNab06282(PDLIM1 Antibody) at dilution of 1:600