

PDPK1 antibody

Product Information

Catalog No.:	FNab06288
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	3-phosphoinositide dependent protein kinase-1
Synonyms:	3-phosphoinositide-dependent protein kinase 1 (hPDK1) PDPK1 PDK1
Observed MW:	48 kDa, 65 kDa
Uniprot ID :	O15530

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200
Image:	



Immunohistochemistry of paraffin-embedded human kidney using FNab06288(PDPK1 Antibody) at dilution of 1:100



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab06288(PDPK1 antibody) at dilution of 1:1000