

## PDXK antibody

### Product Information

Catalog No.:	FNab06294
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Required for synthesis of pyridoxal-5-phosphate from vitamin B6.

### Immunogen information


Immunogen:	pyridoxal(pyridoxine, vitamin B6) kinase
Synonyms:	Pyridoxal kinase Pyridoxine kinase PDXK C21orf124 C21orf97 PKH PNK
Observed MW:	35 kDa, 40 kDa
Uniprot ID :	O00764

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human colon cancer using FNab06294(PDXK antibody) at dilution of 1:50



human liver tissue were subjected to SDS PAGE followed by western blot with FNab06294(PDXK antibody) at dilution of 1:500