

PDZD2 antibody

Product Information

| | |
|---------------|--|
| Catalog No.: | FNab06295 |
| Size: | 100µg |
| Form: | liquid |
| Purification: | Immunogen affinity purified |
| Purity: | ≥95% as determined by SDS-PAGE |
| Host: | Rabbit |
| Clonality: | polyclonal |
| Clone ID: | None |
| IsoType: | IgG |
| Storage: | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.) |

Background


None

Immunogen information

| | |
|--------------|--|
| Immunogen: | PDZ domain containing 2 |
| Synonyms: | PDZ domain-containing protein 2 Activated in prostate cancer protein PDZ domain-containing protein 3 Processed PDZ domain-containing protein 2 PDZD2 AIPC KIAA0300 PDZK3 |
| Observed MW: | None |
| Uniprot ID : | O15018 |

Application

| | |
|-----------------------|---------------------------------|
| Reactivity: | Human, Mouse, Rat |
| Tested Application: | ELISA, IHC, IF |
| Recommended dilution: | IHC: 1:20-1:200; IF: 1:20-1:200 |
| Image: | |



Immunohistochemistry of paraffin-embedded human lung tissue using FNab06295(PDZD2 antibody) at dilution of 1:50