

PEPD antibody

Product Information

Catalog No.:	FNab06307
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Splits dipeptides with a prolyl or hydroxyprolyl residue in the C-terminal position. Plays an important role in collagen metabolism because the high level of iminoacids in collagen.

Immunogen information


Immunogen:	peptidase D
Synonyms:	Xaa-Pro dipeptidase (X-Pro dipeptidase) Imidodipeptidase Peptidase D Proline dipeptidase (Prolidase) PEPD PRD
Observed MW:	55 kDa
Uniprot ID :	P12955

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human lung cancer using FNab06307(PEPD antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab06307(PEPD antibody) at dilution of 1:600