

EIF2AK3 antibody

Product Information

Catalog No.:	FNab06317
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	eukaryotic translation initiation factor 2-alpha kinase 3
Synonyms:	Eukaryotic translation initiation factor 2-alpha kinase 3 PRKR-like endoplasmic reticulum kinase Pancreatic eIF2-alpha kinase (HsPEK) EIF2AK3 PEK PERK
Observed MW:	140 kDa
Uniprot ID :	Q9NZJ5

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human pancreas slide using FNab06317(PERK Antibody) at dilution of 1:50



K-562 cells were subjected to SDS PAGE followed by western blot with FNab06317(PERK antibody) at dilution of 1:1000