

## PGAP3 antibody

### Product Information

|               |  |
|---------------|--|
| Catalog No.:  | FNab06318  |
| Size:         | 100µg  |
| Form:         | liquid   |
| Purification: | Immunogen affinity purified  |
| Purity:       | ≥95% as determined by SDS-PAGE   |
| Host:         | Rabbit   |
| Clonality:    | polyclonal   |
| Clone ID:     | None   |
| IsoType:      | IgG  |
| Storage:      | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.) |

### Background


None

### Immunogen information

|              |   |
|--------------|---|
| Immunogen:   | per1-like domain containing 1   |
| Synonyms:    | Post-GPI attachment to proteins factor 3 COS16 homolog (hCOS16) Gene coamplified with ERBB2 protein PER1-like domain-containing protein 1 PGAP3 CAB2 PERLD1 |
| Observed MW: | 50 kDa  |
| Uniprot ID : | Q96FM1  |

### Application

|                       |                  |
|-----------------------|------------------|
| Reactivity:           | Human            |
| Tested Application:   | ELISA, WB        |
| Recommended dilution: | WB: 1:200-1:1000 |
| Image:                |                  |



HeLa cells were subjected to SDS PAGE followed by western blot with FNab06318( PERLD1 Antibody) at dilution of 1:300