

PRDX2 antibody

Product Information

Catalog No.:	FNab06320
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	peroxiredoxin 2
Synonyms:	Peroxiredoxin-2 Natural killer cell-enhancing factor B (NKEF-B) PRP Thiol-specific antioxidant protein (TSA) Thioredoxin peroxidase 1 Thioredoxin-dependent peroxide reductase 1 Thioredoxin-dependent peroxiredoxin 2 PRDX2 NKEFB TDPX1
Observed MW:	22 kDa
Uniprot ID :	P32119

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab06320(Peroxiredoxin 2 antibody) at dilution of 1:1000