

PEX16 antibody

Product Information

Catalog No.:	FNab06328
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Required for peroxisome membrane biogenesis. May play a role in early stages of peroxisome assembly. Can recruit other peroxisomal proteins, such as PEX3 and PMP34, to de novo peroxisomes derived from the endoplasmic reticulum(ER). May function as receptor for PEX3.

Immunogen information

Immunogen:	peroxisomal biogenesis factor 16
Synonyms:	Peroxisomal membrane protein PEX16 Peroxin-16 Peroxisomal biogenesis factor 16 PEX16
Observed MW:	38 kDa
Uniprot ID :	Q9Y5Y5

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



human liver tissue were subjected to SDS PAGE followed by western blot with FNab06328(PEX16 antibody) at dilution of 1:500