

PEX5L antibody

Product Information

Catalog No.:	FNab06332
Size:	100μg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Accessory subunit of hyperpolarization-activated cyclic nucleotide-gated(HCN) channels, regulating their cell-surface expression and cyclic nucleotide dependence.

Immunogen information


Immunogen:	peroxisomal biogenesis factor 5-like
Synonyms:	PEX5-related protein PEX2-related protein PEX5-like protein Peroxin-5-related protein Peroxisome biogenesis factor 5-like Tetratricopeptide repeat-containing Rab8b-interacting protein (Pex5Rp, TRIP8b) PEX5L PEX5R PXR2
Observed MW:	68 kDa
Uniprot ID :	Q8IYB4

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:2000; IP: 1:200-1:1000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human brain tissue slide using FNab06332(PEX5L Antibody) at dilution of 1:50



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab06332(PEX5L antibody) at dilution of 1:600