

PEX7 antibody

Product Information

Catalog No.:	FNab06333
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Binds to the N-terminal PTS2-type peroxisomal targeting signal and plays an essential role in peroxisomal protein import.

Immunogen information

Immunogen:	peroxisomal biogenesis factor 7
Synonyms:	Peroxisomal targeting signal 2 receptor (PTS2 receptor) Peroxin-7 PEX7 PTS2R
Observed MW:	34 kDa
Uniprot ID :	O00628

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



MDA-MB-453s cells were subjected to SDS
PAGE followed by western blot with
FNab06333(PEX7 antibody) at dilution of 1:800