

## PFDN2 antibody

### Product Information

Catalog No.:	FNab06335
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Binds specifically to cytosolic chaperonin(c-CPN) and transfers target proteins to it. Binds to nascent polypeptide chain and promotes folding in an environment in which there are many competing pathways for nonnative proteins.

### Immunogen information


Immunogen:	prefoldin subunit 2
Synonyms:	Prefoldin subunit 2 PFDN2 PFD2
Observed MW:	17-20 kDa
Uniprot ID :	Q9UHV9

### Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:1000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human lung cancer using FNaB06335(PFDN2 antibody) at dilution of 1:100



L02 cells were subjected to SDS PAGE followed by western blot with FNaB06335(PFDN2 antibody) at dilution of 1:300