

## PFDN4 antibody

### Product Information

Catalog No.:	FNab06336
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Binds specifically to cytosolic chaperonin(c-CPN) and transfers target proteins to it. Binds to nascent polypeptide chain and promotes folding in an environment in which there are many competing pathways for nonnative proteins.

### Immunogen information

Immunogen:	prefoldin subunit 4
Synonyms:	Prefoldin subunit 4 Protein C-1 PFDN4 PFD4
Observed MW:	15 kDa
Uniprot ID :	Q9NQP4

### Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab06336(PFDN4 antibody) at dilution of 1:300