

## PFKL antibody

### Product Information

Catalog No.:	FNab06340
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Catalyzes the phosphorylation of D-fructose 6-phosphate to fructose 1,6-bisphosphate by ATP, the first committing step of glycolysis.

### Immunogen information


Immunogen:	phosphofructokinase, liver
Synonyms:	ATP-dependent 6-phosphofructokinase, liver type (ATP-PFK, PFK-L) 6-phosphofructokinase type B Phosphofructo-1-kinase isozyme B (PFK-B) Phosphohexokinase PFKL
Observed MW:	85 kDa
Uniprot ID :	P17858

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:200-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human liver cancer using FNab06340(PFKL antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab06340(PFKL antibody) at dilution of 1:500