

PGAM2 antibody

Product Information

Catalog No.:	FNab06347
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Interconversion of 3-and 2-phosphoglycerate with 2,3-bisphosphoglycerate as the primer of the reaction. Can also catalyze the reaction of EC 5.4.2.4(synthase) and EC 3.1.3.13(phosphatase), but with a reduced activity.

Immunogen information

Immunogen:	phosphoglycerate mutase 2(muscle)
Synonyms:	Phosphoglycerate mutase 2 BPG-dependent PGAM 2 Muscle-specific phosphoglycerate mutase Phosphoglycerate mutase isozyme M (PGAM-M) PGAM2 PGAMM
Observed MW:	28-29 kDa
Uniprot ID :	P15259

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with FNab06347(PGAM2 antibody) at dilution of 1:1500