

PGBD1 antibody

Product Information

Catalog No.:	FNab06349
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	piggyBac transposable element derived 1
Synonyms:	PiggyBac transposable element-derived protein 1 Cerebral protein 4 PGBD1
Observed MW:	85-90 kDa
Uniprot ID :	Q96JS3

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:200-1:2000; IHC: 1:20-1:200; IF: 1:10-1:100
Image:	



Immunohistochemistry of paraffin-embedded human kidney tissue slide using FNab06349(PGBD1 Antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab06349(PGBD1 Antibody) at dilution of 1:600