

CPQ antibody

Product Information

Catalog No.:	FNab06352
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	plasma glutamate carboxypeptidase
Synonyms:	Carboxypeptidase Q Lysosomal dipeptidase Plasma glutamate carboxypeptidase CPQ LCH1 PGCP
Observed MW:	65-70 kDa
Uniprot ID :	Q9Y646

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



L02 cells were subjected to SDS PAGE followed by western blot with FNab06352(PGCP antibody) at dilution of 1:600