

## PGLYRP3 antibody

### Product Information

Catalog No.:	FNab06357
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Pattern receptor that binds to murein peptidoglycans(PGN) of Gram-positive bacteria. Has bactericidal activity towards Gram-positive bacteria. May kill Gram-positive bacteria by interfering with peptidoglycan biosynthesis. Binds also to Gram-negative bacteria, and has bacteriostatic activity towards Gram-negative bacteria. Plays a role in innate immunity.

### Immunogen information

Immunogen:	peptidoglycan recognition protein 3
Synonyms:	Peptidoglycan recognition protein 3 Peptidoglycan recognition protein I-alpha (PGLYRP1alpha, PGRP-I-alpha) Peptidoglycan recognition protein intermediate alpha PGLYRP3 PGRP1A
Observed MW:	38 kDa
Uniprot ID :	Q96LB9

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



Apoptosised HeLa cells were subjected to SDS PAGE followed by western blot with FNab06357(PGLYRP3 antibody) at dilution of 1:500