

## PGM1 antibody

### Product Information

Catalog No.:	FNab06360
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	0F6
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


This enzyme participates in both the breakdown and synthesis of glucose.

### Immunogen information


Immunogen:	phosphoglucomutase 1
Synonyms:	Phosphoglucomutase-1 (PGM 1) Glucose phosphomutase 1 PGM1
Observed MW:	61 kDa
Uniprot ID :	P36871

### Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:500; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human colon tissue slide using FNab06360( PGM1 Antibody) at dilution of 1:200



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab06360(PGM1 antibody) at dilution of 1:1000