

## PGRMC1 antibody

### Product Information

Catalog No.:	FNab06365
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background


This gene encodes a putative membrane-associated progesterone steroid receptor. The protein is expressed predominantly in the liver and kidney.

### Immunogen information


Immunogen:	progesterone receptor membrane component 1
Synonyms:	Membrane-associated progesterone receptor component 1 (mPR) Dap1 IZA PGRMC1 HPR6.6 PGRMC
Observed MW:	24 kDa
Uniprot ID :	O00264

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200; IF: 1:50 - 1:100
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer using FNab06365(PGRMC1 antibody) at dilution of 1:100



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab06365(PGRMC1 antibody) at dilution of 1:1000