

## EGLN2 antibody

### Product Information

Catalog No.:	FNab06374
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	egl nine homolog 2(C. elegans)
Synonyms:	Prolyl hydroxylase EGLN2 Egl nine homolog 2 Estrogen-induced tag 6 (EIT-6) HPH-3 Hypoxia-inducible factor prolyl hydroxylase 1 (HIF-PH1, HIF-prolyl hydroxylase 1, HPH-1) Prolyl hydroxylase domain-containing protein 1 (PHD1) EGLN2 EIT6
Observed MW:	48 kDa
Uniprot ID :	Q96KS0

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



mouse testis tissue were subjected to SDS PAGE followed by western blot with FNab06374(PHD1 antibody) at dilution of 1:1000