

## EGLN1 antibody

### Product Information

Catalog No.:	FNab06376
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	egl nine homolog 1(C. elegans)
Synonyms:	Egl nine homolog 1 Hypoxia-inducible factor prolyl hydroxylase 2 (HIF-PH2, HIF-prolyl hydroxylase 2, HPH-2) Prolyl hydroxylase domain-containing protein 2 (PHD2) SM-20 EGLN1 C1orf12
Observed MW:	46 kDa
Uniprot ID :	Q9GZT9

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human heart using FNab06376(EGLN1 antibody) at dilution of 1:100



mouse heart tissue were subjected to SDS PAGE followed by western blot with FNab06376(EGLN1 antibody) at dilution of 1:500