

PHF19 antibody

Product Information

Catalog No.:	FNab06383
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)


Background

Immunogen information

Immunogen:	PHD finger protein 19
Synonyms:	PHD finger protein 19 Polycomb-like protein 3 (hPCL3) PHF19 PCL3
Observed MW:	66 kDa
Uniprot ID :	Q5T6S3

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500 - 1:2000
Image:	



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab06383(PHF19 antibody) at dilution of 1:1000