

PHF6 antibody

Product Information

Catalog No.:	FNab06388
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Transcriptional regulator that associates with ribosomal RNA promoters and suppresses ribosomal RNA(rRNA) transcription.

Immunogen information


Immunogen:	PHD finger protein 6
Synonyms:	PHD finger protein 6 PHD-like zinc finger protein PHF6 CENP-31 KIAA1823
Observed MW:	41 kDa
Uniprot ID :	Q8IWS0

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human ovary tumor using FNab06388(PHF6 antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab06388(PHF6 antibody) at dilution of 1:1000