

PHLDA2 antibody

Product Information

Catalog No.:	FNab06397
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background

Plays a role in regulating placenta growth. May act via its PH domain that competes with other PH domain-containing proteins, thereby preventing their binding to membrane lipids(By similarity).


Immunogen information

Immunogen:	pleckstrin homology-like domain, family A, member 2
Synonyms:	Pleckstrin homology-like domain family A member 2 Beckwith-Wiedemann syndrome chromosomal region 1 candidate gene C protein Imprinted in placenta and liver protein Tumor-suppressing STF cDNA 3 protein Tumor-suppressing subchromosomal transferable fragment candidate gene 3 protein p17-Beckwith-Wiedemann region 1 C (p17-BWR1C) PHLDA2 BWR1C HLDA2 IPL TSSC3
Observed MW:	25 kDa
Uniprot ID :	Q53GA4

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000

Image:



human placenta tissue were subjected to SDS PAGE followed by western blot with FNab06397(PHLDA2 antibody) at dilution of 1:800