

## PHLPP2 antibody

### Product Information

Catalog No.:	FNab06399
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	PH domain and leucine rich repeat protein phosphatase-like
Synonyms:	PH domain leucine-rich repeat-containing protein phosphatase 2 PH domain leucine-rich repeat-containing protein phosphatase-like (PHLPP-like) PHLPP2 KIAA0931 PHLPPL
Observed MW:	147 kDa
Uniprot ID :	Q6ZVD8

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human colon cancer tissue slide using FNab06399(PHLPPL Antibody) at dilution of 1:50



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with FNab06399(PHLPPL Antibody) at dilution of 1:600