

## phospho-Marcks (Ser167/170) antibody

### Product Information

Catalog No.:	FNab06404
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	myristoylated alanine rich protein kinase C substrate
Synonyms:	Myristoylated alanine-rich C-kinase substrate (MARCKS) Protein kinase C substrate, 80 kDa protein, light chain (80K-L protein, PKCSL) MARCKS MACS PRKCSL
Observed MW:	80 kDa
Uniprot ID :	P29966

### Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:3000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab06404( Marcks Antibody) at dilution of 1:500