

PHOX2A antibody

Product Information

Catalog No.:	FNab06407
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May be involved in regulating the specificity of expression of the catecholamine biosynthetic genes. Acts as a transcription activator/factor. Could maintain the noradrenergic phenotype.

Immunogen information

Immunogen:	paired-like homeobox 2a
Synonyms:	Paired mesoderm homeobox protein 2A ARIX1 homeodomain protein Aristaless homeobox protein homolog Paired-like homeobox 2A PHOX2A ARIX PMX2A
Observed MW:	35-40 kDa
Uniprot ID :	O14813

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with FNab06407(PHOX2A Antibody) at dilution of 1:600