

## PHYHIP antibody

### Product Information

Catalog No.:	FNab06413
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Its interaction with PHYH suggests a role in the development of the central system.

### Immunogen information

Immunogen:	phytanoyl-CoA 2-hydroxylase interacting protein
Synonyms:	Phytanoyl-CoA hydroxylase-interacting protein Phytanoyl-CoA hydroxylase-associated protein 1 (PAHX-AP1, PAHXAP1) PHYHIP DYRK1AP3 KIAA0273
Observed MW:	35 kDa, 43-45 kDa
Uniprot ID :	Q92561

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab06413(PHYHIP antibody) at dilution of 1:1000