

PIK3CB antibody

Product Information

Catalog No.:	FNab06419
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	phosphoinositide-3-kinase, catalytic, beta polypeptide
Synonyms:	Phosphatidylinositol 4,5-bisphosphate 3-kinase catalytic subunit beta isoform (PI3-kinase subunit beta, PI3K-beta, PI3Kbeta, PtdIns-3-kinase subunit beta) Phosphatidylinositol 4,5-bisphosphate 3-kinase 110 kDa catalytic subunit beta (PtdIns-3-kinase subunit p110-beta, p110beta) Serine/threonine protein kinase PIK3CB PIK3CB PIK3C1
Observed MW:	110 kDa
Uniprot ID :	P42338

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IF, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:10-1:100
Image:	



Immunohistochemistry of paraffin-embedded human lung cancer tissue slide using FNab06419(PIK3CB Antibody) at dilution of 1:200 heat mediated antigen retrieved with Tris-EDTA buffer(pH9).



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab06419(PIK3CB antibody) at dilution of 1:500