

## PIGO antibody

### Product Information

Catalog No.:	FNab06443
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Ethanolamine phosphate transferase involved in glycosylphosphatidylinositol-anchor biosynthesis. Transfers ethanolamine phosphate to the GPI third mannose which links the GPI-anchor to the C-terminus of the proteins by an amide bond(By similarity).

### Immunogen information


Immunogen:	phosphatidylinositol glycan anchor biosynthesis, class O
Synonyms:	GPI ethanolamine phosphate transferase 3 Phosphatidylinositol-glycan biosynthesis class O protein (PIG-O) PIGO
Observed MW:	74 kDa
Uniprot ID :	Q8TEQ8

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human liver using FNab06443(PIGO antibody) at dilution of 1:100



HeLa cells were subjected to SDS PAGE followed by western blot with FNab06443(PIGO antibody) at dilution of 1:300