

## PNN antibody

### Product Information

Catalog No.:	FNab06459
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	pinin, desmosome associated protein
Synonyms:	Pinin 140 kDa nuclear and cell adhesion-related phosphoprotein Desmosome-associated protein Domain-rich serine protein (DRS protein, DRSP) Melanoma metastasis clone A protein Nuclear protein SDK3 SR-like protein PNN DRS MEMA
Observed MW:	82 kDa, 140 kDa
Uniprot ID :	Q9H307

### Application

Reactivity:	Human
Tested Application:	ELISA, IHC, IF, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human kidney using FNab06459(PNN antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab06459(PNN Antibody) at dilution of 1:600