

PINX1 antibody

Product Information

Catalog No.:	FNab06460
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Microtubule-binding protein essential for faithful chromosome segregation. Mediates TRF1 and TERT accumulation in nucleolus and enhances TRF1 binding to telomeres. Inhibits telomerase activity. May inhibit cell proliferation and act as tumor suppressor.

Immunogen information


Immunogen:	PIN2-interacting protein 1
Synonyms:	PIN2/TERF1-interacting telomerase inhibitor 1 Liver-related putative tumor suppressor Pin2-interacting protein X1 Protein 67-11-3 TRF1-interacting protein 1 PINX1 LPTL LPTS
Observed MW:	40 kDa
Uniprot ID :	Q96BK5

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human skin tissue slide using FNab06460(PINX1 Antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab06460(PINX1 antibody) at dilution of 1:1000